Supplementary Table 1. Primers used for segregation testing and *in vitro* splice assay

Name	Primer sequence (5'-3')	Amplicon size (bp)
hu PALM3 Ex4-5 XhoI F	aattctcgagTCTGGGACAGTGATGCTGTT	864 bp
hu PALM3 Ex4-5 BamHI R	attggatccTGACTTGCCAAGACCACAGA	
OTOA Ex19 F	GAATGTCTTTGGGATGAGCTCTTGG	378 bp
OTOA Ex19 R	TTTAGAGTAAGCTTCAAGACCCACC	
TMC1 Ex 16-17 F	CCTAGTTTCTCCCTTGTGCCT	643 bp
TMC1 Ex16-17 R	AATTCAGAGCCAGCACACAG	
SD6 F	TCTGAGTCACCTGGACAACC	NA
SA2 R	ATCTCAGTGGTATTTGTGAGC	